

# FEVAR is the Best Treatment for Juxtarenal Aneurysms But There are Limitations



**Thomas L. Forbes, MD**

*R. Fraser Elliott Chair in Vascular Surgery*

*Peter Munk Cardiac Centre, University Health Network*

*University of Toronto*

# DISCLOSURE

## Thomas L. Forbes, MD

- **Research Grants:** Cook Medical – Alpha Abdominal Graft



# FENESTRATED STENT GRAFTS






Figure 3

# MISALIGNMENT

**Vertical**





**Horizontal**




# MISALIGNMENT

## Stent Graft Planning




### NON STANDARD DEVICE REQUEST

PAGE 2 of 2









# MISALIGNMENT - ROTATION



## What are the short term effects of intraoperative fenestration misalignment in with FEVAR?


1. *What is the incidence of vertical and horizontal misalignment?*
2. *Is misalignment associated with poor clinical outcomes?*

# HORIZONTAL MISALIGNMENT





# VERTICAL MISALIGNMENT



# RESULTS

**Table 1.** Patient Demographics and Comorbidities [Mean  $\pm$  SEM or No (%)]

	Aligned [N=29]	Misaligned [N=36]	P-Value
<b>Demographics</b>			
Age	74 $\pm$ 2	78 $\pm$ 1	0.02
Male Gender	27 (93)	26 (72)	0.05
Body Mass Index	29.8 $\pm$ 1.1	29.4 $\pm$ 1.2	
<b>Comorbidities</b>			
Hypertension	24 (82)	30 (83)	> 0.5
Dyslipidemia	22 (76)	30 (83)	> 0.5
Chronic Renal Failure	11 (37)	12 (33)	> 0.5
Coronary Artery Disease	16 (55)	25 (69)	0.3
Stroke	5 (17)	0 (0)	0.01
MI or Prior PCI/CABG	13 (45)	16 (44)	> 0.5
COPD	13 (45)	16 (44)	> 0.5
Type 2 Diabetes	10 (34)	8 (22)	0.4
Peripheral Vascular Disease	7 (24)	15 (42)	0.18
Tobacco use	25 (86)	28 (77)	> 0.5
Pervious Abdominal Surgery	5 (17)	4 (11)	> 0.5
Previous Aortic Surgery	2 (7)	0 (0)	0.19

# RESULTS

**Table 4.** Procedural details and 30-Day/In-hospital post-operative clinical outcomes

	Aligned [N=29]	Misaligned [N=36]	P-Value
<b>Total Procedural Time (min)</b>	262 ± 15	351 ± 25	0.007
<b>Mean Fluoroscopy Time (min)</b>	85 ± 6	106 ± 8	0.04
<b>Mean Contrast Volume (ml)</b>	193 ± 13	239 ± 20	0.07
<b>Hospital Length of Stay (days)</b>	5 ± 0.5	10 ± 1.5	0.005
<b>30-day Complications</b>			
Atrial fibrillation	1 (3)	2 (6)	> 0.5
Paraplegia	1 (3)	4 (11)	0.37
Myocardial infarction	2 (7)	2 (6)	> 0.5
Ischemic Colitis	0 (0)	4 (11)	0.12
Renal Failure	0 (0)	1 (3)	> 0.5
Death	0 (0)	4 (11)	0.12
<b>No. of Endoleaks</b>			
Type 1	0 (0)	0 (0)	> 0.5
Type 2	4 (14)	2 (6)	0.39
Type 3	3 (10)	2 (6)	> 0.5
<b>Intraoperative Target Vessel Cannulation Failure</b>	0 (0)	4 (11)	0.12
<b>Intraoperative Target Vessel Complication <sup>1</sup></b>	1 (3)	8 (8)	0.04
<b>End-organ Ischemia and/or Death</b>	1 (3)	11 (31)	0.008

# MISALIGNMENT

- >50% of patients had some degree of horizontal or vertical misalignment intraoperatively
- Patients with misalignment had
  - Higher incidence of composite adverse events (31% vs 3%)
  - Longer operative and fluoroscopy times
  - Higher rates of intraoperative target vessel complications

A non-planning cause of misalignment is rotation

Are there predictors of FEVAR Rotation?

1. What are the short-term clinical outcomes in patients with intra-operative stent graft rotation?
2. Are there anatomical markers of the arterial geometry that can predict stent graft rotation?

# ILIAC ARTERIES


## – Local geometric properties

- **Calcification**
- **Radius**
  - Maximum Inscribed Sphere Radius
- **Curvature**

$$\kappa(s) = \frac{|c'(s) \times c''(s)|}{|c'(s)|^3}$$

- **Torsion**

$$\tau(s) = \frac{[c'(s) \times c''(s)] \cdot c'''(s)}{|c'(s) \times c''(s)|^2}$$



# GEOMETRIC VARIABLES

## TORTUOSITY

$$\chi = L/D$$



# GEOMETRIC VARIABLES

## TORSION



Measures the local deviation of the centerline from the osculating plane





Measures how sharply the centerline twists out of the plane of curvature

$$\tau = \frac{\det(r', r'', r''')}{\|r' \times r''\|^2} = \frac{(r' \times r'') \cdot r'''}{\|r' \times r''\|^2}$$

Torus-Knoten mit Tangentialvektor (braun), Normalenvektor (grün) und Binormalenvektor (blau)







$$R = 6.03 \left( \sum Torsion (mm^{-1}) \right) + 0.0276 (Calcium Volume (mm^3)) + 0.221 (DeviceLength (mm)) - 101.4$$

# FEVAR ROTATION

- Patients who experience stent graft rotation have significantly higher rates of severe post-operative complications (36% vs 0%)
- FEVAR rotation causes fenestration misalignment and is predictable
- Iliac artery torsion and calcification are associated with stent graft rotation



## WHAT CAN A SURGEON DO TO MINIMIZE DEVICE ROTATION?

1. To develop a mechanically realistic aortoiliac bench top model
2. To validate quantitative anatomical markers which have been observed clinically
3. To evaluate the effect of operator technique on stent graft rotation

- Flexible models
  - Molds are 3D printed in PLA at a resolution of 0.2 mm
  - Polyvinyl alcohol cryogel is then cast using a 15% solution with 4 freeze thaw cycles
  - Effective rigidity is controlled by altering the thickness of the model walls
- Idealized models
  - Helically structured to have a constant torsion
  - Torsion specified by altering the pitch and radius of the helix
- Patient specific models
  - Segmented from pre-operative CT angiograms



## EXPERIMENTAL APPARATUS







Patient 1



Patient 2




Patient 3



**Table 4-1.** Measured rotation during ZFEN deployment in patient specific aortoiliac phantoms compared with the observed intraoperative rotation during patient deployment.

Patient No.	Model Rotation (°)	Intraoperative Rotation (°)	Absolute Error (°)	Iliac Torsion (mm <sup>-1</sup> )	Iliac Calcium Volume (mm <sup>3</sup> )
1	2.3	0.2	2.1	5.2	413
2	27	29	2	13.1	1352
3	34	30	4	8.2	1244



**Figure 4-4.** Evaluation of the effect of operator insertion technique in rigid idealized models at varying levels of torsion. **Straight:** No correction of device orientation during insertion. **Multiple:** Insert the device, note orientation, fully remove the device, correct the orientation accordingly, and re-insert. **Corrected:** Gradually correct the orientation of the device as needed during insertion.



# FEVAR MISALIGNMENT & ROTATION

- Torsion and rigidity/calcification combined are direct causes of intraoperative stent graft rotation
- In-vivo correction of orientation significantly increases the observed rotation
  - If clinically safe, fully remove the device, adjust the orientation and reinsert the device



## Supervisors

Cristina Amon, ScD  
Thomas Forbes, MD

## Students

Sean Crawford  
Ryan Sanford  
Patrick Lin  
Fabian Alvarado

Helen Genis  
Gabrielle Sebaldt  
Michael Lancaster  
Masha Itkina  
Wanyu Zhang

## Collaborators

Graham Roche-Nagle (UHN)  
Naomi Eisenberg (UHN)  
Elrasheed Osman (UHN)  
Leonard Tse (UHN)  
Glen Van Arsdell (SickKids)  
Osami Honjo (SickKids)  
Lucy Roche (UHN)  
Richard Weisel (UHN)  
Haruki Ide (SickKids)  
Christoph Haller (SickKids)  
Gilles Soulez (CHUM)  
Claude Kauffmann (CHUM)  
David Roy (CHUM)  
Hossein Mohammadi (CHUM)  
Elena Di Martino (Calgary)

**TED ROGERS CENTRE  
FOR HEART RESEARCH**



UNIVERSITY OF TORONTO  
FACULTY OF APPLIED SCIENCE & ENGINEERING

